



ANDOVER HEALTHCARE



Animal Health

A tradition of excellence in product innovation.

TABLE OF CONTENTS

COMPANY PROFILE..... 3

PRODUCT FEATURES..... 4

PETFLEX..... 6-7

PETFLEX NL..... 8

COFLEX VET..... 9

POWERFLEX..... 10-11

PAWTAPE..... 12

PET FLEX AFD..... 14

POWERFLEX AFD..... 15



COMPANY PROFILE

Andover Healthcare is a leading manufacturer of high quality cohesive bandages and tapes for the animal health, healthcare and sports medicine industries. Founded in 1976, Andover is committed to new and innovative technology and distributes its patented products around the world. All Andover products are made in the USA.

Andover strives to develop break-through products to meet the ever-changing needs of the animal health industry. Product development is Andover's specialty and customer service is our commitment. At Andover Healthcare, it's not just about our products, technology or service - it's about all of these things, all of the time.

Veterinary professionals, equine trainers and pet owners alike have embraced the patented technology of PetFlex and PowerFlex for the support it provides for their animals.



Product Features

FLEXIBLE COHESIVE BANDAGES

Andover offers a variety of cohesive flexible bandage options.

- ✓ Available latex or latex-free
- ✓ Self-adhering - sticks to itself, not hair or skin
- ✓ Easy Hand-Tear patented technology - no scissors needed
- ✓ Controlled compression - will not constrict over time
- ✓ Available in a variety of sizes and strengths



COHESIVE BANDAGES

	Latex-Free	Easy Hand-Tear Application	Patented Technology	Variety of Colors	Fun Prints	Bitter Taste No Chew	Small Animal Use	Large Animal Use
Coflex-Vet				●			●	●
PetFlex		●	●	●	●	●	●	
PetFlex_{NL}	●	●	●	●			●	
PowerFlex		●	●	●	●	●		●

WOUND CARE

ADHESIVE TAPE

PetFlex •AFD	Absorbent foam dressing for small animal use	PawTape [™]	Synthetic adhesive veterinary tape
PowerFlex •AFD	Absorbent foam dressing for large animal use		



Cohesives + Tapes



PetFlex®

A FLEXIBLE COHESIVE BANDAGE

EasyTear
PATENTED TECHNOLOGY



Contains
Natural
Rubber
Latex

SMALL ANIMAL



Give your pet
protection, support,
& *sparkle* with:
PetFlex®



A flexible, cohesive bandage with quick
and easy application

EasyTear
PATENTED TECHNOLOGY

U.S. Patent No. 5,762,623
European Patent No. 0839222

Key Features:

- ✓ Quick easy application - no scissors necessary
- ✓ EasyTear® patented technology
- ✓ Controlled compression - will not constrict
- ✓ Strong - 15 lbs. tensile strength
- ✓ Sweat & water resistant - stays in place
- ✓ Available in colors, fun prints, glitter & bitter-tasting No Chew

Uses for
PetFlex®:

- ✓ Light Compression
- ✓ Fixing
- ✓ Wound Care

PetFlex® is available in single
color cases or in color packs
for variety!



When ordering, please indicate catalog # and color code. For example, 2100TN-030 (XX = color code) Most products are available in bulk packaging, call for details.

Catalog #	Size	Rolls/Case	Colors/Prints	Catalog #	Size	Rolls/Case	Colors/Prints
2100XX-030	1" x 5 yds.	30		2300XX-024	3" x 5 yds.	24	
2200XX-036	2" x 5 yds.	36		2400XX-018	4" x 5 yds.	18	

PetFlex[®] NL



A FLEXIBLE COHESIVE BANDAGE



SMALL ANIMAL



Give your pet non-latex protection & support with:
PetFlex[®] NL



A non-latex flexible, cohesive bandage with quick and easy application



U.S. Patent No. 5,762,623, 6,156,424
European Patent No. 0839222

Color Pack (CP)



Includes: Neon Pink, Blue, Purple, Light Blue, Neon Green, Red

Key Features:

- ✓ Strong - 12 lbs. tensile strength
- ✓ EasyTear[®] patented technology
- ✓ Quick easy application - no scissors necessary
- ✓ Sweat & water resistant - stays in place
- ✓ Controlled compression - will not constrict
- ✓ Available in a variety of colors

PetFlex[®] NL is available in single color cases or in a color pack for variety!



Uses for PetFlex^{NL}: ✓ Light Compression ✓ Fixing ✓ Wound Care

When ordering, please indicate catalog # and color code. For example, 6200TN-036 (XX = color code) Most products are available in bulk packaging, call for details.

Catalog #	Size	Rolls/Case	Colors/Prints	Catalog #	Size	Rolls/Case	Colors/Prints
6100XX-030	1" x 5 yds.	30	●○●●●●●●	6300XX-024	3" x 5 yds.	24	●○●●●●●●
6200XX-036	2" x 5 yds.	36	●○●●●●●●	6400XX-018	4" x 5 yds.	18	●○●●●●●●



CoFlex® Vet

A FLEXIBLE COHESIVE BANDAGE

LARGE + SMALL



Economical protection & identification with: CoFlex® Vet



A flexible, cohesive bandage with quick and easy application

Key Features:

- ✓ Economical self-adhering bandage
- ✓ Provides protection & superior adhesion, easy removal
- ✓ Easy to use - apply quickly in a continuous motion
- ✓ Available in neon colors to help identify infected cows
- ✓ Works well even in cold weather
- ✓ Won't constrict - no repeat visits to remove constricting bandages
- ✓ Sweat & water resistant - stays in place

Uses for CoFlex Vet: ✓ Fixing Bandage ✓ Light Compression ✓ Wound Care ✓ Identification

Color Pack (CP) Rainbowpack (RB)



CoFlex Vet is available in single color cases or in color packs for variety!

Teal (TE)	Includes: Neon Pink, Blue, Purple, Light Blue, Neon Green, Red	Includes: Teal, White, Blue, Purple, Red, Green
Black (BK)	Light Blue (LB)	Orange (OR)
Purple (PU)	Yellow (YL)	Red (RD)
Neon Green (NG)	Green (GR)	White (WH)
Neon Pink (NP)	Blue (BL)	Tan (TN)

When ordering, please indicate catalog # and color code. For example, 3510TN-030 (XX = color code) Most products are available in bulk packaging, call for details.

Catalog #	Size	Rolls/Case	Colors/Prints	Catalog #	Size	Rolls/Case	Colors/Prints
3510XX-030	1" x 5 yds.	30		3540XX-018	4" x 5 yds.	18	
3520XX-036	2" x 5 yds.	36		3540XX-100	4" x 5 yds.	100	
3530XX-024	3" x 5 yds.	24		3560XX-012	6" x 5 yds.	12	

PowerFlex®

A FLEXIBLE COHESIVE BANDAGE



Contains Natural Rubber Latex

LARGE ANIMAL



Protection, support, & sparkle for equines: PowerFlex®



A strong, flexible cohesive bandage ideal for support. High abrasion resistance keeps it from shredding.

Key Features:

- ✓ Strong - 23 lbs. tensile strength
- ✓ Won't shred - high abrasion resistance
- ✓ Straight hand-tear
- ✓ Quick, easy application - no scissors needed
- ✓ Controlled compression - will not constrict
- ✓ Sweat & water resistant - won't saturate or slip
- ✓ Available in colors, fun prints, glitter & No Chew

Uses for PowerFlex®: ✓ Ligament Support ✓ Light Compression
✓ Fixing Bandage ✓ Wound Care

PowerFlex® is available in single color cases or in color packs for variety!

U.S. Patent No. 5.762.623
European Patent No. 0839222

Purple Glitter (GU)	Blue Glitter (GB)
Pink Glitter (GN)	Black Glitter (GK)
Maroon (MR)	Red Glitter (GD)
Grey (GY)	Teal (TE)
Black (BK)	Yellow (YL)
Light Blue (LB)	Green (GR)
Orange (OR)	Blue (BL)
Purple (PU)	Red (RD)
Neon Green (NG)	White (WH)
Neon Pink (NP)	Tan (TN)

When ordering, please indicate catalog # and color code. For example, 3840TN-018 (XX = color code) Most products are available in bulk packaging, call for details.

Catalog #	Size	Rolls/Case	Colors/Prints	Catalog #	Size	Rolls/Case	Colors/Prints
3840XX-018	4" x 5 yds.	18		3840XX-100	4" x 5 yds.	100	

PawTape™

VETERINARY ADHESIVE TAPE



Contains
Natural
Rubber
Latex

SMALL ANIMAL



Veterinary adhesive tape with *fun* paw prints:
PawTape™



A synthetic adhesive veterinary tape that provides protection + consistent support.

Key Features:

- ✓ Strong - 33 lbs. tensile strength
- ✓ Synthetic adhesive tape
- ✓ Crisp, easy hand-tear
- ✓ Delivers consistent support
- ✓ Easy to apply
- ✓ Made in the USA

Uses for PawTape: ✓ Immobilization ✓ Splinting + Strapping
 ✓ Patient Positioning ✓ Hold Dressings in Place



Catalog #	Size	Rolls/Case	Colors/Prints	Catalog #	Size	Rolls/Case	Colors/Prints
291-005-288	½" x 10 yds.	288	●	291-020-072	2" x 10 yds.	72	●
291-010-144	1" x 10 yds.	144	●				



Wound Care



PetFlex® AFD

ABSORBENT FOAM DRESSING

LATEX
Contains
Natural
Rubber
Latex

SMALL ANIMAL



All-in-One wound management system: **PetFlex® AFD**

 All-in-one foam pad dressing & bandage

Key Features:

- ✓ Quick and easy to apply
- ✓ Absorbs the equivalent of 29 gauze pads (4" size)
- ✓ Controlled compression - will not constrict
- ✓ EasyTear® Technology - no scissors needed
- ✓ Conforms to wound bed - reduces maceration and promotes healing
- ✓ Won't stick to wound bed

Uses for PetFlex^{AFD}: ✓ Puncture Wounds ✓ Abrasions ✓ Lacerations

EasyTear
PATENTED TECHNOLOGY

U.S. Patent No. 5,762,623
European Patent No. 0839222

Applying PetFlex AFD (Leg Wound)



After cleaning the wound, place foam pad over wound.

Place foam pad over wound.







Cover wound with foam pad & unroll to expose absorbent foam pad. Place foam pad over wound. Cover wound with foam pad & wrap bandage around so it covers first layer. Prior to each wrap around leg, pull bandage out approx. 6".



Continue to wrap at desired tension with 50% overlap.



Seal end of bandage to bottom layer with finger pressure for a strong, solid seal.

Catalog #	Size	Rolls/Case	Colors/Prints	Catalog #	Size	Rolls/Case	Colors/Prints
2220PW-016	2" x 2.5 yds.	16		2440PW-008	4" x 2.5 yds.	8	
2220RPW-016 (Reverse Wind)	2" x 2.5 yds.	16		2440RPW-008 (Reverse Wind)	4" x 2.5 yds.	8	

PowerFlex® AFD

ABSORBENT FOAM DRESSING

LATEX
Contains
Natural
Rubber
Latex

LARGE ANIMAL



All-in-One wound management system: PowerFlex® AFD

 All-in-one foam pad dressing & bandage

Key Features:

- ✓ Foam pad and PowerFlex cohesive in one
- ✓ Quick and easy to apply
- ✓ Absorbs and holds up to 10 times its weight in fluid
- ✓ Strong, flexible PowerFlex base
- ✓ Controlled compression - will not constrict
- ✓ Straight, hand-tear - no scissors needed
- ✓ Effective on puncture wounds, abrasions and lacerations
- ✓ Conforms to wound bed - reduces maceration and promotes healing
- ✓ Won't stick to wound bed

Applying PowerFlex AFD (Leg Wound)

U.S. Patent No. 5.762.623
European Patent No. 0839222



After cleaning the wound, unroll to expose absorbent foam pad.

Place foam pad over wound.



Cover wound with foam pad & wrap bandage around so it covers first layer. Prior to each wrap around leg, pull bandage out approx. 6".



Continue to wrap at desired tension with 50% overlap.



Seal end of bandage to bottom layer with finger pressure for a strong, solid seal.

Uses for PowerFlex_{AFD} ✓ Puncture Wounds ✓ Abrasions ✓ Lacerations

Catalog #	Size	Rolls/Case	Colors/Prints
3744BK-008	4" x 2.5 yds.	8	●
3744RBK-008 (Reverse Wind)			

Catalog #	Size	Rolls/Case	Colors/Prints
3744S-008 (Sterile)	4" x 2.5 yds.	8	○



Scan code to watch a video on how to apply PowerFlex





ANDOVER HEALTHCARE, INC.
9 FANARAS DRIVE, SALISBURY, MA 01952
PHONE: 978.465.0044 TOLL FREE: 800.432.6686
FAX: 978.462.0003

www.andoverhealthcare.com